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**DialInspect-OSM** is an automated controlled microscope, especially designed for size, shape and segment analysis of abrasives and tools. It covers the whole range of particle sizes from micron up to coarse grit.

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## Hardware:

Main part of the system is the professional microscope with 5 magnification steps, digitally controlled lighting installation and high-quality digital industrial video camera. The microscope is equipped with a motor driven Z-axis and motor driven mechanical staging equipment for X- and Y-axis. All axes are controlled either by joystick or directly by the control program. The video camera is connected via USB3 to an image acquisition system, which makes the pictures available to the image processing software.

## Software:

The systems operation is controlled by an easy to handle program interface, which runs on a usual PC.

## Principle of operation:

A number of particles has to be spread onto a glass plate. When inspecting particles which are few micron in size wet application of the powder is recommended. The glass plate is moved in scan motion by the XY-stage. The array size is automatically chosen depending on the actual magnification level. The image acquisition part acquires and analyses consecutive pictures from every part of the glass plate, so a very large number of single particles can be inspected with a high image resolution.

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Grain size : 1µm...3mm

Scan range : 76 x 52 mm (XY-stage)

Processing time : app. 15 minutes  
(preparation and test for 2000 particles, size 50µm)

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## Calculated parameters:

More than 20 parameters were calculated per particle.

- Minimum feret diameter
- Maximum feret diameter
- Circle equivalent grain diameter
- Total area
- Perimeter
- Convex perimeter
- Moments of inertia
- Feret elongation
- Compactness
- Roughness
- Ellipticity
- Percentage of the "transparent" appearing area
- Relative brightness of the "transparent" appearing area
- Colour coordinates in the CIE L\*a\*b-system
- ...

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## Report :

The control program generates printed reports which contain statistical results for every of the parameters along with histograms and density plots. For scientific work on the acquired data a transfer function copies every number for all of the particles into an EXCEL workbook.

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## Special features:

The control program displays a data grid, where one row contains the parameter of one single particle. A simple click onto any row let the picture of the described particle appear in a display. So the correlation of a particles picture and the numbers which describe the appearance numerically becomes easier to understand.

A special add-on is the virtual sorting facility. It allows to apply a virtual sorting process on the inspected particles: sorting ranges can be specified by random arrays of two of the above parameters. A density plot shows the distribution of the particles along the chosen parameters. In which box the program has to put the particles from within an area of the density plot you can specify simply by dragging the mouse. After "sorting" the program shows you the amount of particles in every box and builds up one picture for every box, which shows the particles now inside.